


Disposable Blood Pressure Sensor**Low Cost Disposable****Dielectric Gel Barrier****Fully Tested & Calibrated****mV or 24 bits I2C Output****Descriptions**

The SA5650 is a fully piezo resistive silicon pressure sensor for use in medical applications such as filtration, chromatography or invasive blood pressure monitoring. The sensor is designed to be used with automated assembly equipment and can be dropped directly into a customer's pressure housing. The sensor is designed to meet the requirements as described in the Association for the Advancement of Medical Instrumentation (AAMI) specification BP22 for medical application like Blood Pressure Transducers.

The pressure sensor consists of a pressure sensing element mounted on a PCB board. SMT resistors on the PCB are laser trimmed or ASIC adjacent to the sensing element for compensation and calibration. A polycarbonate or polysulfone plastic cap is attached to the PCB substrate to provide an easy method of attachment to the customer's assembly and protection for the sensing element. A dielectric gel is placed over the sensor to provide electrical and fluid isolation.

The products are shipped on a tape and reel. Performance characteristics and packaging can be easily tailored on a special order basis to meet the requirements of specific customers.

Features and Applications

Low Cost

Small Size and Reliable Performance

Gel Isolation for Liquids

Operates from 5°C to 50°C

Compatible for Automated Assembly

1% Accuracy for Replacements

5.0 or 40 $\mu\text{V/V/mmHg}$ Sensitivity

Customization for OEM Applications

Surgical Procedures

Intensive Care Units

Infusion Pumps

Kidney Dialysis Machines

Vacuum Assisted Birth

Intrauterine Monitoring

Filtration

Chromatography

Performance specifications

Unless otherwise specified: Supply Voltage: 6.0 V_{DC} , Ambient Temperature: 23°C

Operating Pressure Range	-50		300	mmHg	
	1		150	psi	
Over Pressure	125		150	psi	
Zero Pressure Offset	-20		20	mmHg	
Sensitivity	4.95	5.00	5.05	$\mu\text{V/V/mmHg}$	1
Calibration	97.5	100	102.5	mmHg	2
Linearity and Hysteresis (-50 to 100 mmHg)	-1		1	mmHg	3
Linearity and Hysteresis (100 to 300 mmHg)	-1.5		1.5	% Output	3
Linearity and Hysteresis (5 to 150 psi)	-1.5		1.5	% Output	3
Input Impedance					
SA5650-mV	1200		3200	Ω	
SA5650-DO	NA		NA	Ω	
Output Impedance	300		1000	Ω	


Output Symmetry	0.95		1.05	Ratio	4
Supply Voltage	1	6	10	V _{DC} or V _{AC rms}	
Risk Current (@ 230 or 120 V _{AC rms} , 50 or 60Hz)			2	µA	
Warm-Up Time		5		Seconds	
Frequency Response		1200		Hz	
Offset Drift			1	mmHg	5
Thermal Span Shift	-0.1		0.1	%/°C	6
Thermal Offset Shift	-0.3		0.3	mmHg/°C	6
Phase Shift (@ 5KHz)			5	Degrees	
Light Sensitivity (3000 Foot Candle)		1		mmHg	
Defibrillator withstand (400 joules)	5			Discharges	7
Sterilization (ETO)	5			Cycles	8
Operating Temperature	5		50	°C	
Storage Temperature	-25		+70	°C	
Humidity (External)	10-90% (non-condensing)				
Operating Product Life	168 Hours for Medical related application				
Shelf Life	5 Years				
Dielectric Breakdown	10000V _{DC}				
Media Interface	Dielectric Gel				
Volume Displacement	4.5 x 10 ⁻⁴ mm ³				
Weight	0.5 grams				

Notes


1. Sensitivity can be calibrated according to customer request from 5-40uV/mmHg
2. Output of sensor with no pressure applied and a certain specified resistor placed across +SUPPLY to +OUTPUT.
3. Best fit straight line.
4. Defined as common mode symmetry between any output and supply terminal.
5. Over an 8 hour time period after a 10 second warm-up.
6. Over operating temperature range 5°C–50°C with respect to 23°C.
7. One discharge per minute performed by customer.
8. Sterilization performed by customer, compatible with ETO, GAMMA, or E-Beam sterilization.
9. After curing, meniscus of gel shall be flush to .035" [.89mm] below surface
10. The configuration with polycarbonate plastic cap is not suitable for surface mount operation. The configuration with polysulfone plastic cap is suitable for surface mount operation.

Equivalent Circuit

EQUIVALENT CIRCUIT





SA5650-mV



SA5650-DO


Dimensions and Pad configurations

Units: Inch (millimeter)



DETAIL B
SCALE: 1:1
PADS CONFIGURATION

SA5650-mV




DETAIL B
SCALE: 1:1
PADS CONFIGURATION


SA5650-DO

Shipping Tape and Reel information

UNITS: INCHES



HEAT PRESSURE STICKING





- 1) THE PACKAGE CAN BE CUSTOMIZED ACCORDING TO CUSTOMER REQUIREMENT.
- 2) EACH REEL PACKED WITH MAX 1440 SENSORS, HAVE A DIFFERENT LOT NO.
- 3) NO BLANK SPACES ON THE TAPE BETWEEN SENSORS.
- 4) REELS PACKED IN CARTON BOXES. EACH BOX CONTAIN 6 REELS. PACKED IN 2 PE BAGS.
- 5) INNER BAG SEALED WITH TRANSPARENT TAPE.
- 6) EACH REEL LABELED WITH P.N, LOT NO. AND QTY INCLUDED.
- 7) CARTON BOX AND INNER BAG LABELED WITH NAME OF SA, SA P.N. PRODUCT DESCRIPTION, ORDER NO. QTY IN EACH REEL.

Ordering information

SA5650	-	XX	-	XXX	-	X
		Output type- mV		Pressure range		Pressure Unit
				030		M- mmHg
		Output type- I2C		100		P- PSI
				150		

Digital Output Calculation



I2C Protocol

I2C INTERFACE PARAMETERS					
PARAMETERS	SYMBOL	MIN	TYP	MAX	UNITS
SCLK CLOCK FREQUENCY	F _{SCL}	100		400	KHz
START CONDITION HOLD TIME RELATIVE TO SCL EDGE	t _{HDSTA}	0.1			μs
MINIMUM SCL CLOCK LOW WIDTH @1	t _{LOW}	0.6			μs
MINIMUM SCL CLOCK HIGH WIDTH @1	t _{HIGH}	0.6			μs
START CONDITION SETUP TIME RELATIVE TO SCL EDGE	t _{SUSTA}	0.1			μs
DATA HOLD TIME ON SDA RELATIVE TO SCL EDGE	t _{HDDAT}	0			μs
DATA SETUP TIME ON SDA RELATIVE TO SCL EDGE	t _{SUDAT}	0.1			μs
STOP CONDITION SETUP TIME ON SCL	t _{SUSTO}	0.1			μs
BUS FREE TIME BETWEEN STOP AND START CONDITION	t _{BUS}	2			μs

@1 COMBINED LOW AND HIGH WIDTHS MUST EQUAL OR EXCEED MINIMUM SCL PERIOD.

I2C INTERFACE TIMING DIAGRAM

